

COMPLETE MOLECULAR SOLUTION FOR SARS-COV-2



For the detection of SARS-CoV-2 by PCR

! All products are CE marked

DETECTION OF SARS-COV-2, FLU & RSV

Triplex	✓ SARS-CoV-2: orf1ab and N gene
	✓ Flu A/B
	✓ RSV A/B
	✓ Endogenous control with RNase-P
	✓ Only 1 MasterMix
	✓ New to the market

Product code:

v2 ABCOW6 (100 reactions)
v4 ABCOWI (100 reactions),
Triplex ABSFR3 (100 reactions)
Multiplex ABSCR3 (100 reactions)

DETECTION OF SARS-COV-2 ONLY

SARS-CoV-2 kit v2	✓ orf1ab (including RdRp)
	✓ E gene
	✓ Internal control
	✓ Extensively used in UK and France
SARS-CoV-2 kit v4	✓ orf1ab (including RdRp)
	✓ E gene
	✓ N gene
	✓ Endogenous control with RNase-P
	✓ New to the market

DETECTION OF SARS-COV-2 + MULTIPLE RESPIRATORY PATHOGENS (15 TARGETS)

SARS-CoV-2 respiratory panel	✓ SARS-CoV-2: orf1ab and N-gene	✓ Metapneumovirus
	✓ Flu A/B	✓ Rhinovirus
	✓ RSV A/B	✓ <i>Legionella pneumophila</i>
	✓ Adenovirus	✓ <i>Mycoplasma pneumoniae</i>
	✓ Enterovirus	✓ Parainfluenza 1 - 4
	✓ Endogenous control with RNase-P	
	✓ Only 3 MasterMixes	
	✓ New to the market	

ABACUS

Diagnostica

Sample to answer instrument, the GenomEra, to detect SARS-CoV-2

! All products are CE marked

Product code:
CDX-130-01-40



Dual target: RdRp and E genes

No RNA extraction needed

Results in 70min

A SARS-CoV-2/Flu/RSV assay is in development

vircell

MICROBIOLOGISTS

A third party control
for SARS-CoV-2 assays

AmpliRun® Total SARS-CoV-2 control

! All products are CE marked

Product code:
MBTC030-R



External quality control for nucleic acids detection

Validation of the complete analytical process:
extraction & amplification

Complete inactivated microorganism in a matrix



For further information please call
+44 (0) 1474 874 426
or visit our website

launchdiagnostics.com